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UKH 602

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Circulating cooler

UKH 602

Operating Instructions

FRYKA
Kältetechnik

Ohmstraße 4
D- 73730 Esslingen
Fon: +49 (0)711 310599-0
Fax: +49 (0)711 310599-29
E-Mail: info@fryka.de
Web: www.fryka.de

Inhalt // Index

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

en // Operating Instructions Circulating Cooler UKH 602

Compliance with the operating instructions is a basic prerequisite for operating the FRYKA cooling devices safely and for achieving the specified product qualities and service features. FRYKA-Kältetechnik GmbH is not liable for any bodily injuries, material or financial damages due to non-compliance with the operating instructions. There is no liability for material defects in such cases.

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1 General information

1.1 Intended Use

This equipment is designed for temperature control and transportation of certain liquids. An external circuit is made at the connections. The pressure for the external circuit must be determined and may not be exceeded.

DANGER: DANGER OF DEATH DUE TO EXPLOSION

- › Do not install the device in a hazardous area.
- › Do not take the device in a hazardous area in operation.

WARNING

- › Inappropriate use can result in considerable personal injury or material damages. As a inappropriate use is, among other things: direct temperature control of foods, or the using as a medical device and operation without coolant or an unsuitable coolant.

1.2 Scope of supply

- ◆ Circulating Cooler
- ◆ Power cord
- ◆ Operating instructions
- ◆ Hose connector (10mm) for drain cock
- ◆ Four hose connectors (2x 10mm/ 2x 8mm) for coolant connections

1.3 Explanation of the safety notes

Safety notes are marked by a pictogram and a signal word. The signal word describes the severity of the threatened danger.

DANGER

- › Imminent threat to the life and health of persons (severe injury or death).

WARNING

- › Potential threat to the life and health of persons (severe injury or death).

CAUTION

- › Potentially dangerous situation (slight injuries or material damages).

NOTE

- › Obligation to behave or act in a certain way in handling the unit.

1.4 Conformity

The refrigeration cycle has been checked for leakage. The unit complies with the relevant safety regulations as well as EC directives 2006/95/EC, 2004/108/EC and 2011/65/EU.

2 Safety and warning notes

2.1 Before you start operation of the unit

Read these operating instructions carefully! They contain important information about installation, safe operation and maintenance of the circulating cooler. Keep all documents for later use. Make certain that those responsible for the unit and users have read the operating instructions completely and understood them.

Third parties should not make any changes to the unit. If any changes are made without the manufacturer's agreement, all explanations about the circulating cooler lose their validity.

WARNING

- ▶ Improper interventions or repairs may lead to considerable personal injury or material damages.
- ▶ Only skilled personnel trained by the manufacturer or authorised technicians may make changes or repairs or conduct maintenance.

The manufacturer does not assume any liability for damages due to technical changes to the circulating cooler, improper handling, abuse or use of the circulating cooler without heeding the operating instructions.

2.2 General safety instructions

WARNING

- ▶ This unit is not designed for persons with physical, sensory or mental impairments, or persons without enough experience or knowledge, unless they are shown how to use the circulating cooler and are initially supervised by a person responsible for their safety.
- ▶ Pull out the power cord for maintenance, cleaning or repair work. Always pull on the plug, not on the cable.
- ▶ Remove a defective unit from the power supply immediately and never use any defective unit. Check the unit and accessories for damage at regular intervals.

NOTE

- ▶ Never cover or close the ventilation grille of the circulating cooler.

3 Installation and set-up

3.1 Install the circulating cooler

DANGER: FIRE AND ELECTRIC SHOCK RISKS DUE TO MOISTURE

- › If current-carrying parts or the power cord become moist, a short circuit may occur. The unit is designed for use in enclosed rooms. Do not operate the unit outside or in areas moist or splashed with water.

The circulating cooler is suitable for installation only in dry and ventilated/air-conditioned rooms. Heed the permissible *ambient temperature* (Chapter "*Technical data*"). The circulating cooler should not be left in direct sunshine or be near a source of heat such as radiators.

CAUTION: COLLISION/SQUASHING DUE TO ROLLING OR FALLING OF THE DEVICE BY MISHANDLING

- › Place the device with two people on a flat and straight ground.
- › Make sure that the device is in a safe and stable position.
- › Do not tilt the unit.
- › Pay attention to a non-slip area with sufficient load-capacity.
- › Lock if necessary the brakes at the rolls.

NOTE

- › Do not operate the unit in a tilted position. The refrigeration unit may otherwise fail to work.
- › Too much heat at the installation site may have the same effect.

Install the unit in such a way as to ensure optimum ventilation. Maintain a distance of at least 30 cm from the ventilation grille to the nearest object (wall, etc.).

3.2 Connect the circulating cooler

Wait at least one hour after connecting the circulating cooler before starting it up. Check the unit for damages each time before starting the circulating cooler up.

DANGER: FIRE/ELECTRIC SHOCK DUE TO FAULTY CONNECTION

- › Only connect the circulating cooler to a properly installed socket with earthing contact.
- › The socket may be earthed with a maximum 16A slow-to-blow earth.
- › Never operate a device that has damaged parts or a damaged mains connection.
- › Damaged cables can be exchanged by a qualified electrician.

NOTE

- › Operate the circulating cooler only with the mains voltage, current and frequency indicated on the rating plate. Not heeding the instructions may cause the unit to malfunction or cause defects.

WARNING: FIRE- AND OVERHEATING HAZARD

- › Do not use any extension cords or distributor banks. A risk of fire and overheating.

3.3 Connect coolant hoses

CAUTION: INJURY FROM USING UNSUITABLE/DEFECTIVE TUBES/HOSE CONNECTION

- › Use only for the range of temperatures and pressures suitable coolant thermally insulated hoses and hose connections.
- › Check the coolant hoses in regular intervals for material fatigue and replace them if necessary.
- › Prevent kinking of a hose.
- › Safeguard the connections against slipping by using hose clips.

WARNING. DANGER OF BURSTING

- › Determine before commissioning the max. admissible pressure for the external circuit. This is determined by the weakest link (for example glass apparatuses). Make sure that this is not exceeded.

Now set up the connection with your external system by using the coolant hoses. To do this, connect the cooling unit of your system to the entry (upper plug), then from the exit (medium plug) back to your system. Check the tightness and the tightness of the hose clamps.

3.4 Fill cooling unit and first turn

DANGER: RISK OF FIRE / ELECTRIC SHOCK HAZARD BY HUMIDITY IN THE DEVICE

- › Do not fill the device over a water pipe – the pressure destroys the heat exchanger.
- › Do not pour coolant on the device.

CAUTION:

- › Follow the instructions of the safety data sheet of the cooling medium if necessary, in particular the instructions for use of personal protective equipment (PPE).
- › Do not move the device in the filled state, in order to avoid spillage of the coolant.

Fill the cooling system through the filler plug with a funnel and an appropriate coolant. Make sure that the coolant is suitable for the temperature range (freezing). Turn on the pump during the filling process, as the system must be bled before placing into operation. Pump and chiller go only above a certain minimum level in operation. The filled coolant is now pumped into the circuit externally connected.

NOTE:

- › Observe the notes for **coolant** (Chapter "**specifications**"). Non-observance will damage the device.
- › Never let the pump run without coolant.

After filling turn on the device and set the desired temperature on the controller.

4 Control and operation

4.1 Control

Press the main switch to turn on the circulating cooler. The current interior temperature is shown in the display. The desired interior temperature is displayed by pushing the "set" key.

To change the desired interior temperature press and hold the "set" key and use the arrow-keys to set the desired temperature. On completion of the adjustment process, check the setpoint by pressing the SET button.

The LED „E1“ lights up permanently when the minimum level of the coolant tank has been reached and no coolant alarm is present.

4.2 Alarm

Status and error messages are shown in the display. An acoustic signal also indicates an alarm. To deactivate the acoustic signal push the arrow-down-key.

You can find the precise significance of error codes in the *Status and error messages table* (Chapter "*Status and error messages*").

4.3 Switching off

Turn off the circulating cooler by pushing the main switch.

5 Design and function

The refrigeration unit of the cooler consists of a hermetically sealed refrigerant circuit, which is ventilated with an axial fan. The heat exchanger is made of non-ferrous-metal and has a 3/4" filler plug and two 12mm hose connections with a 10mm inner thread. There is also a drain plug (lowest level connection), through which one the system can be drained completely. In the cooling device a circulation pump is installed, which is not self priming, but can be operated for filling and sealing of the heat exchanger with open systems, as long as the open bath is higher than the cooling device. The system can be vented at the vent screw on the filler plug. In this case, the system should be filled by the external bath.

6 Maintenance

6.1 Refrigeration unit

Clean the fin coil heat exchanger of the cooling device every six months whatever the degree of soiling at the installation site. The fin coil heat exchanger is located on the back side of the unit.

NOTE

- › A soiled fin coil heat exchanger results in overheating and outage of the unit!

CAUTION

- › Avoid direct contact with the fin coil heat exchanger – the fin coils are sharp-edged. Danger of cuts!

Other components of the cooling device are designed for continuous operation and do not require any maintenance.

7 Troubleshooting

WARNING

- ▶ Only technicians authorised by the manufacturer should do repairs or intervene in the refrigeration cycle or the control of the device.
- ▶ Improper interventions or repairs can result in considerable bodily injury and/or material damages.
- ▶ Pull out the mains power plug before performing maintenance or repair work.

7.1 Internal thermal protection

The compressors of the refrigeration unit are equipped with internal thermal protection. This switches off the compressors and consequently the cooling if there is overheating. As a result, the interior temperature increases and a temperature alarm follows.

Check the following possible causes:

- ◆ Is the fin coil heat exchanger soiled?
- ◆ Is ambient temperature too high?
- ◆ Is optimal ventilation of the cooling device ensured?
- ◆ Are the ventilators operating?

The cooling is switched on again automatically after the compressor has cooled down (approx. one to two hours).

You can obtain an overview of the status and error messages in the following table.

You can obtain additional information and help when there are malfunctions and error messages from the manufacturer.

7.2 Block/unblock setpoint adjustment

The keyboard lock enables blocking the controller keys. When locked it is not possible to change the desired temperature. Please contact the manufacturer to activate or deactivate the key lock in the parameters.

7.3 Status and error messages

Error display	Cause	Remedy
F1L	Sensor failure: short circuit	Sensor defect; send the unit to the manufacturer for repair.
F1H	Sensor failure: sensor break	Sensor defect; send the unit to the manufacturer for repair.
F3L	Boundary alarm: temperature to low	Electronic circuit defect; send the unit to the manufacturer for repair.
F3H	Boundary alarm: temperature to high, cooling device overheated	<ul style="list-style-type: none"> - Clean (Chapter "Refrigeration unit") the fin coil heat exchanger. - Ambient temperature too high (>30° C)? - Is optimal ventilation of the cooling device ensured?
---	Keyboard lock active	See Block/unblock setpoint <i>adjustment</i> (Chapter " Block/unblock setpoint adjustment ")
Flashing display	Temperature alarm: the current temperature lies outside the working range of the cooling device.	<ul style="list-style-type: none"> - Clean (Chapter "Refrigeration unit") the fin coil heat exchanger. - Check the ventilator. - Ambient temperature too high (>30° C)? - Is optimal ventilation of the cooling device ensured? - Switch off the cooling device, let it stand for 1-2 h and then switch it on again
EP	Data loss in EE-Prom	Controler defect; send the unit to the manufacturer for repair.

8 Decommissioning, disposal

8.1 Decommissioning

Turn off the unit as described.

Leave the mains power cord disconnected if the cooling device is not used for a longer time (energy saving).

8.2 Disposing of old cooling devices

Old circulating coolers are not valueless waste. Valuable raw materials can be recycled with environmentally compatible disposal. Have a specialist firm dispose of your old circulating cooler, or contact your supplier or manufacturer. The circulating cooler is to be disposed of in waste separate from unsorted residential waste (domestic refuse).

Do not damage the refrigeration cycle of an old circulating cooler during transport, to avoid uncontrolled leakage of the refrigerant and oil still in it.

9 Transport, packaging and storage

CAUTION

- › Bring the device to room temperature before emptying.
- › Move and transport the device only completely emptied.

NOTE

- › Transport the circulating cooler only in an upright position and well packed to prevent damages.
- › Only ship the circulating cooler on pallets via a forwarding company.

9.1 Disposing of packaging

The packaging protects your circulating cooler against transport damages. Please be helpful: Dispose of the packaging in an environmentally compatible way. Your municipal authorities can provide information about current disposal options.

10 Technical data

Operation	-10°C to +40°C
Control accuracy	+/- 1 K
Cooling capacity in watts	
at +20°C	570
at +10°C	500
at -10°C	60
Circulating pump	
Output	12l/min
Discharge head	6mWs
External dimension [WxDxH] [cm]	31 x 36 x 52
Weight [kg]	32
Ambient temperature	+12°C to +30°C
Electrical connection	230V / 50Hz
Max. current [A]	2,8
Coolant tank	2,6l

10.1 Coolant

The used refrigerant must be compatible with non-ferrous metals (copper/brass), to avoid damages on the heat exchanger. Do not use distilled or deionized / demineralised water. Due to its corrosive nature, it leads to corrosion in the cooling water circuit. Because of these properties, demineralized water or seawater is also unsuitable. Do not use iron-containing water or water with a high lime content (hard water), this leads to rusting or calcification in the coolant circuit. As not suitable coolant apply cooling water with suspended solids and untreated, not purified water.


NOTICE:

- › The use of an improper coolant will damage the device.

10.2 Parts list

Code	Description	Article number
M1	Compressor	530FR11GX230
M2	Condenser unit	611FCEV100D
M3	Pump 230V / 50 Hz	705MLYE660550P
W1	Condenser	611FCEV100D
W2	Heat exchanger	630WTDLK602
Y1	Capillary	-
Ft	Filter dryer	660FT1562
S1	Switch Cooling	790S820A230
S2	Switch Pump	790S820A230
S3	Float switch	730KSS
R1	Controller	710ST71IGBA
Ks	Low heat device socket with fuses	-
F 1	Temperature sensor controller	790PT10006MJ
K1	Release relay	740MR4W230
Contains refrigerant R134a (GWP 1300). See type plate for filling capacities.		

10.3 Wiring diagramm

Klemmleiste:

PE	N	L	N1	L1	N2	L2	L3	N4	L4	x2	x3	x4	x5	x6	PE

10.4 Refrigeration circuit

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© FRYKA-Kältetechnik GmbH. All rights reserved.