



















# Hei-MIX Platform Shakers – multifarious possibilities

SHAKERS							ACCESSORIES							
DEVICE			Speed (rpm)	Orbit/stroke (mm)	Usable area (mm)	Loading capacity max. (kg)								
							Perforated platform	Erlenmeyer clamps <sup>1</sup> max. configuration	Erlenmeyer flask attachment <sup>1,2</sup>	Separatory funnel clamps <sup>1,2</sup> max. configuration	Attachment for separatory funnel <sup>1</sup>	Tension roller attachment	Multi-tier attachment	
							pc.	ml	ml	pc.	ml	ml		
ORBITAL	Rotamax 120		20 – 300	20	220 x 220	2	100	16 8 5 3	25/50 100 250 500	–	–	–	100	–
	Unimax 1010		30 – 500	10	290 x 258	5	1000	20 14 8 4 2	25/50 100 250 500/1000 2000	25 – 1000	–	50 – 100	1000	–
	Unimax 2010		20 – 400	20	390 x 340	10	2000	36 23 12 9 5 3	25/50 100 250 500 1000 2000	–	4 3 3 2	250 500 1000 2000	–	2000
VIBRATING	Vibramax 100		150 – 1.350	3	220 x 220	2	100	16 8 5 3	25/50 100 250 500	–	–	–	100	–
	Titramax 100/101		150 – 1.350	1,5/3	220 x 220	2	–	–	–	–	–	–	–	–
	Titramax 1000		150 – 1.350	1,5	290 x 258	5	–	–	–	–	–	–	–	–
WAVE	Polymax 1040		2 – 50	–	290 x 258	5	1000	20 14 8 4 2	25/50 100 250 500/1000 2000	25 – 1000	–	50 – 100	1000	–
	Polymax 2040		2.5 – 50	20	390 x 340	10	2000	36 23 12 9 5 3	25/50 100 250 500 1000 2000	–	4 3 2	250 500/1000 2000	–	2000
ROCKING	Duomax 1030		2 – 50	–	290 x 258	5	1000	20 14 8 4 2	25/50 100 250 500/1000 2000	25 – 1000	–	50 – 100	1000	–
RECIPROCATING	Promax 1020		30 – 250	32	290 x 258	5	1000	20 14 8 4 2	25/50 100 250 500/1000 2000	25 – 1000	–	50 – 100	1000	–
	Promax 2020		20 – 400	20	390 x 340	10	2000	36 23 12 9 5 3	25/50 100 250 500/1000 2000	–	4 3 2	250 500/1000 2000	–	2000

# Hei-MIX Platform Shakers – multifarious possibilities

SHAKER				ACCESSORIES		
OVERHEAD	DEVICE	Speed (rpm)	Loading capacity max. (kg)	Adapter		
	Reax 2 	20 – 100	1	 Rack for 20 Falcon tubes/test tubes á 15 ml		
OVERHEAD	DEVICE	Speed (rpm)	Loading capacity max. (kg)	Adapter (max. number)		
	Reax 20 	1 – 16 2 – 32	30	Bottle height mm	pc.	ml
				290 – 250	–	–
				248 – 208	1	1000
				206 – 166	2	1000
				202 – 162	1	500
				160 – 120	1	1000
				160 – 120	1	500
				114 – 74	2	500
				72 – 32	1	1000
				72 – 32	1	500