

Intelligent Evaporation

MADE IN
GERMANY

The average operational **lifespan of 10 years** is backed by a **3 year warranty** and makes your purchase a worthwhile investment.

Hei-VAP Industrial + Hei-VOLUME Distimatic Pro = Hei-VAP³ Pro Industrial Package Saves time – increases throughput – reduces costs

Expandable with various peripherals included in the Workstation AddOn Industrial package.

Cost-effective alternative to 50 l and 100 l evaporator

A factor of 3 times higher throughput compared to an 8 h working day in manual operation

Automatic filling and refilling of the rotary flask

Automatic distillate drainage through overpressure

Automatic drainage of the residue through overpressure (only the 24 / 7 packages)

Automatic shutdown of all system components at the end of the process and in the event of a disturbance

Automatic cleaning program for glass components and media-carrying hoses

Complete efficiency and security in unattended and continuous use

Calibration-free sensors reliably detect any medium

Enabling and maintaining constant conditions for higher throughput

Safe processing of high volumes – even unattended day / night

Usable evaporator system and full flexibility of configuration

Guided menu navigation of the Hei-VOLUME Distimatic Pro automatic module for smooth setting and storage of all process parameters

Multilingual user interface

Individual process setting through selectable operating modes (sensor, time, manual)

Automatic backfeed of the heating bath

Easy cleaning of the heating bath thanks to concealed heating coil and drain valve

Space-saving solution with fixed place for all system components on the rotary evaporator

Base cart (optional) enables easy accessibility to the device from all sides

