

PRECISE HEATING & STIRRING HEI-PLATE MIX 'N' HEAT MAGNETIC STIRRER

HIGH-RESISTANT HEATING PLATE

REFRACTORY HOUSING

800 W

MADE IN
GERMANY

3-year warranty on
all devices and an op-
erational lifespan of
more than 10 years

Hei-PLATE Mix 'n' Heat Ultimate

Interface for more precision and comfort

For the Pt1000 or the double sensor Pt1000 to measure the temperature exactly until inside the liquid

RS 232, LAN, USB A and USB B interfaces for datalogging as well as the remote control and monitoring via PC

Handling of the electric Heidolph lab jack

Well-designed features for a carefree working

Different heating modes (FAST, PRECISE, PRECISE+) and safety features handle your samples optimally

Three-stage setting of the stirring speed (start-up behavior). Additionally options for the rotational direction and the stir bar detector

Ramp programming and creating profiles for a precise and quick reproducibility of your results

Intuitive handling for a quick start

Standby button in view to switch the device easily on and off

Display and timer per main function (stirring and heating) to have an eye on all values and settings at anytime

The locking function of the push-and-rotary controls with LED rings protects against accidental changing of the values

User-friendly menu navigation for the activation/deactivation of well-designed features
Language versions: German, English



ACCESSORIES RECOMMENDATION



Lab jack



Double temperature sensor



Heating and bath attachments



High-performance condenser