




America

CERTIFICATE

No. U8 069850 0012 Rev. 00

Holder of Certificate: **Heidolph Instruments GmbH & Co. KG**
Walpersdorfer Str. 12
91126 Schwabach
GERMANY

Certification Mark:



C US

Product: **Laboratory Equipment**
(Benchtop Series)

This product was voluntarily tested to the relevant safety requirements referenced on this certificate. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL and a Standards Council of Canada accredited Certification body.

Test report no.: 713204515

Date, 2021-07-16

(Benedikt Pulver)



America

CERTIFICATE

No. U8 069850 0012 Rev. 00

Model(s): Hei-VAP Value, Hei-VAP Value Digital, Hei-VAP Advantage, Hei-VAP Precision, Hei-VAP Core, Hei-VAP Expert, Hei-VAP Ultimate
 Heizbad Hei-VAP (Heating bin)

Tested according to: UL 61010-1:2012/R:2016-04
 CAN/CSA-C22.2 No. 61010-1:2012/U2:2016-04
 UL 61010-2-010:2015
 CAN/CSA-C22.2 No. 61010-2-010:2015

Model Name	Hei-VAP Value, Hei-VAP Value Digital, Hei-VAP Advantage, Hei-VAP Precision, Hei-VAP Core, Hei-VAP Expert, Hei-VAP Ultimate	Heizbad Hei-VAP (Heating bin)
Rated voltage:	115 or 230 VAC	115 or 230 VAC
Rated frequency:	50/60Hz	50/60Hz
Rated power:	Max. 1400 W	Max. 1300 W
Protection class:	I	I