

Viton® extension tube for Cassette small - wt: 0.9 | id: 0.89 | od: 2.69 (per meter)

P/N: 036303730



Extension Tubing for Cassette small

- Application: Excellent acid resistance at high temperatures
- Features:
 - Low gas permeability
 - Resistant to solvents and corrosives at high temperatures
- Material: Fluorocarbon rubber, thermoformed Viton® B (67% fluorinated); opaque black
- Complies with the standards: GLP conform
- Temperature range: -30 to +205 °C
- Sterilization: Sterilization is not recommended
- Restriction: Limited performance life

Used with:

Acids:

excellent

Premium Laboratory Equipment

Lyes:	excellent
Solvents:	varies; tests recommended
Pressure:	good
Vacuum:	good
Viscous media:	good
Sterile media:	satisfactory

Viton® extension tube for Cassette small - wt: 0.9 | id: 0.89 | od: 2.69 (per meter) - Technical Data

Inner diameter	0.89 mm
Outer diameter	2.69 mm
Wall thickness (wt)	0.9 mm
Max. Pressure (continuous / short time)	0.5 / 1.5 bar
Suction height	7 mH ₂ O
SP quick Advantage 06 / Precision 06 (min/max)	not suitable
SP quick Value 06 (min/max)	not suitable
SP quick Advantage 01 / Precision 01 (min/max)	not suitable
SP quick Value 01 (min/max)	not suitable
SP standard/vario Advantage 06 / Precision 06 (min/max)	not suitable
SP standard/vario Value 06 (min/max)	not suitable
SP standard/vario Advantage 01 / Precision 01 (min/max)	not suitable
SP standard/vario Value 01 (min/max)	not suitable
Hei-FLOW Advantage 01 / Precision 01 - pump head C 12 - Cassette small (min/max)	0.03 – 1 ml/min
Hei-FLOW Advantage 01 / Precision 01 - pump head C 4 - Cassette small (min/max)	0.24 – 6 ml/min
Hei-FLOW Value 01 - pump head C 12 - Cassette small (min/max)	0.10 – 1 ml/min
Hei-FLOW Value 01 - pump head C 4 - Cassette small (min/max)	0.57 – 6 ml/min
Hei-FLOW Value 01 - pump head C 8 - Cassette medium (min/max)	not suitable
Hei-FLOW Value 01 - pump head C 8 - Cassette large (min/max)	not suitable
Hei-FLOW Advantage 01 / Precision 01 - pump	not suitable

head C 8 - Cassette medium (min/max)

Hei-FLOW Advantage 01 / Precision 01 - pump not suitable

head C 8 - Cassette large (min/max)

Flow rates pertain to

Tygon® (Standard) tubing and water