

Viton® extension tube for Cassette small - wt: 0.9 | id: 2.79 | od: 4.59 (per meter)

P/N: 036303750



Extension Tubing for Cassette small

- Application: Excellent acid resistance at high temperatures
- Features:
 - Low gas permeability
 - Resistant to solvents and corrosives at high temperatures
- Material: Fluorocarbon rubber, thermoformed Viton® B (67% fluorinated); opaque black
- Complies with the standards: GLP conform
- Temperature range: -30 to +205 °C
- Sterilization: Sterilization is not recommended
- Restriction: Limited performance life

Used with:

Acids:	excellent
--------	-----------

Premium Laboratory Equipment

Lyes:	excellent
Solvents:	varies; tests recommended
Pressure:	good
Vacuum:	good
Viscous media:	good
Sterile media:	satisfactory

Viton® extension tube for Cassette small - wt: 0.9 | id: 2.79 | od: 4.59 (per meter) - Technical Data

Inner diameter	2.79 mm
Outer diameter	4.59 mm
Wall thickness (wt)	0.9 mm
Max. Pressure (continuous / short time)	0.5 / 1.5 bar
Suction height	7 mH ₂ O
SP quick Advantage 06 / Precision 06 (min/max)	not suitable
SP quick Value 06 (min/max)	not suitable
SP quick Advantage 01 / Precision 01 (min/max)	not suitable
SP quick Value 01 (min/max)	not suitable
SP standard/vario Advantage 06 / Precision 06 (min/max)	not suitable
SP standard/vario Value 06 (min/max)	not suitable
SP standard/vario Advantage 01 / Precision 01 (min/max)	not suitable
SP standard/vario Value 01 (min/max)	not suitable
Hei-FLOW Advantage 01 / Precision 01 - pump head C 12 - Cassette small (min/max)	0.29 – 9 ml/min
Hei-FLOW Advantage 01 / Precision 01 - pump head C 4 - Cassette small (min/max)	2 – 36 ml/min
Hei-FLOW Value 01 - pump head C 12 - Cassette small (min/max)	0.69 – 8 ml/min
Hei-FLOW Value 01 - pump head C 4 - Cassette small (min/max)	4 – 37 ml/min
Hei-FLOW Value 01 - pump head C 8 - Cassette medium (min/max)	not suitable
Hei-FLOW Value 01 - pump head C 8 - Cassette large (min/max)	not suitable
Hei-FLOW Advantage 01 / Precision 01 - pump	not suitable

head C 8 - Cassette medium (min/max)

Hei-FLOW Advantage 01 / Precision 01 - pump	not suitable
head C 8 - Cassette large (min/max)	

Flow rates pertain to	Tygon® (Standard) tubing and water
-----------------------	------------------------------------