

PharMed® two-stop tubing for Cassette small - wt: 0.9 | id: 0.51 | od: 2.31

P/N: 036303850



Two-stop tubing (40 cm) for Cassette small

- Application: Ideal for medical, lab and research uses
- Features:
 - High fatigue strength under repeated reversed bending stresses
 - Non-toxic, biocompatible
 - Very low gas permeability
 - Well suited to acids and bases
- Material: Thermoplastic elastomer on a polypropylene base with plasticizers; excellent tensile strength; opaque
- Complies with the standards: USP Class VI, GLP, Pharmacopoea and Europaea
- Temperature range: -51 to +135 °C
- Sterilization: Can be autoclaved or sterilized with ethylene oxide or sterilized by irradiation
- Restriction: Additives may migrate

Used with:

Premium Laboratory Equipment

| | |
|----------------|------------|
| Acids: | good |
| Lyes: | good |
| Solvents: | unsuitable |
| Pressure: | good |
| Vacuum: | excellent |
| Viscous media: | good |
| Sterile media: | excellent |

PharMed® two-stop tubing for Cassette small - wt: 0.9 | id: 0.51 | od: 2.31 - Technical Data

| | |
|--|---------------------|
| Inner diameter | 0.51 mm |
| Outer diameter | 2.31 mm |
| Wall thickness (wt) | 0.9 mm |
| Max. Pressure (continuous / short time) | 0.5 / 1,5 bar |
| Suction height | 7 mH ₂ O |
| SP quick Advantage 06 / Precision 06 (min/max) | not suitable |
| SP quick Value 06 (min/max) | not suitable |
| SP quick Advantage 01 / Precision 01 (min/max) | not suitable |
| SP quick Value 01 (min/max) | not suitable |
| SP standard/vario Advantage 06 / Precision 06 (min/max) | not suitable |
| SP standard/vario Value 06 (min/max) | not suitable |
| SP standard/vario Advantage 01 / Precision 01 (min/max) | not suitable |
| SP standard/vario Value 01 (min/max) | not suitable |
| Hei-FLOW Advantage 01 / Precision 01 - pump head C 12 - Cassette small (min/max) | 0.01 – 0.54 ml/min |
| Hei-FLOW Advantage 01 / Precision 01 - pump head C 4 - Cassette small (min/max) | 0.08 – 2 ml/min |
| Hei-FLOW Value 01 - pump head C 12 - Cassette small (min/max) | 0.02 – 0.42 ml/min |
| Hei-FLOW Value 01 - pump head C 4 - Cassette small (min/max) | 0.17 – 2 ml/min |
| Hei-FLOW Value 01 - pump head C 8 - Cassette medium (min/max) | not suitable |
| Hei-FLOW Value 01 - pump head C 8 - Cassette large (min/max) | not suitable |
| Hei-FLOW Advantage 01 / Precision 01 - pump | not suitable |

head C 8 - Cassette medium (min/max)

| | |
|---|--------------|
| Hei-FLOW Advantage 01 / Precision 01 - pump | not suitable |
| head C 8 - Cassette large (min/max) | |

| | |
|-----------------------|------------------------------------|
| Flow rates pertain to | Tygon® (Standard) tubing and water |
|-----------------------|------------------------------------|