

Hei-VAC Vario Station

P/N: 036308310



Fully controllable stand-alone pumping unit including vacuum controller.
For Hei-VAP Value and Hei-VAP Advantage

- Precise vacuum control avoids bumping and eliminates foaming of your evaporation solution
- Due to automatic vacuum supply to the process parameters you have more time for important laboratory activities

Condenser is not in scope of delivery.

Hei-VAC Vario Station - Technical Data

Weight	6 kg
Diaphragm pump	Three-stage
Suction capacity	1.7 m³/h
Ultimate vacuum	5 mbar
Combinable with a condenser	yes
Power input	160
Dimensions without condenser (w/d/h)	294 x 192 x 299 mm