

BR 14 Collapsible-Blade Impeller

P/N: 036300370



Especially when using narrow-necked vessels, e.g. round-bottom flasks, an stirring tool with collapsible blades is used. Just like our collapsible-blade impeller BR 14.

This blade impeller is ideal for dissolving and homogenizing in vessels with the following characteristics:

- ideal vessel diameter: 180 mm
- ideal ground clearance: 9 mm
- ideal filling level: 180 mm
- ideal volume: 4,600 ml

The collapsible-blade impeller BR 14 has a total length of 400 mm and is made of stainless steel. Its stirring shaft has a diameter of 8 mm.

BR 14 Collapsible-Blade Impeller - Technical Data

Blade size	90 x 10 mm
Material	stainless steel V4A
Length	400 mm
Shaft diameter	8 mm
Max. rotation speed	800 rpm