

## TR 21 Radial-Flow Impeller

P/N: 036300390



### Radial-Flow Impeller

- These impellers are recommended for applications which require average speed

- For mixing tasks with little or average viscosity up to  $< 500 \text{ mPa s}$
- Ideal for gassing of liquids
- These impellers create a radial flow

## TR 21 Radial-Flow Impeller - Technical Data

Blade size	50 mm
Material	stainless steel V4A
Length	400 mm
Shaft diameter	8 mm
Max. rotation speed	2,000 rpm